

Procedural skills training

Powered by

**TOTUM**<sup>®</sup>

Accessible | Affordable | Functional



# Totum Digitising the learning of procedural skills

With the Totum Mobile App, direct procedural skills training can now be recorded, shared, and assessed from any location. Totum mobile allows the digitization of any procedural skills trainer, allowing programmes to fast track their digitalisation of procedural skills teaching and assessment.

The user's mobile phone is mounted above task trainers such as the Inovus Basic Surgical Skills Base. The phone records procedures, which can then be reviewed and assessed in the Totum web application.

## Remote learning and certification



# Powered by Totum for purposeful practice

Totum's video capture and debrief feature allows users to review their performance and receive remote coaching from connected mentors. Mentors and examiners can issue certificates of competence directly through the Totum portal.

By linking with the Totum platform, the app allows over 30 direct procedural skills to be assessed remotely by trainers. In this way coaching can be closely focused on the individual's technique and specific areas for improvement.

The Totum Mobile App can be used to capture training sessions on any brand of procedural skills trainer by setting the mobile phone up in its cradle and focusing the phone's camera on the skills trainer. Procedural videos are captured and aggregated in Totum, allowing programs to easily and cost-effectively digitise their procedural skills learning, and assessment.

## Key features of Totum

- Designed for use with iOS and Android devices
- Recorded procedure simulations stored, managed, and available to share remotely
- Procedures categorised by skill type, making them easily retrievable by mentors and peers for review and assessment
- Create certificate templates specific to your programme and issue these within Totum
- Expand the capacity of your skills lab by delivering remote skills training and assessment

[Totum mobile licence](#)

Product code:  
00001589



**TOTUM** To learn more about the complete Totum ecosystem of connected surgical training, visit [www.inovus.org](http://www.inovus.org)



## Surgical Training Remastered

---

Global Headquarters  
St Helens, United Kingdom  
E: [info@inovus.org](mailto:info@inovus.org)  
T: +44 (0)1744 752 952  
www: [inovus.org](http://inovus.org)

USA Headquarters  
Tampa, Florida USA  
E: [info@inovus.org](mailto:info@inovus.org)  
T: 512 481 7691  
www: [inovus.org](http://inovus.org)

