


bozzini®

Hysteroscopy Simulator

Accessible | Affordable | Functional



Bozzini

Hysteroscopy Simulator

The Bozzini Hysteroscopy Simulator is a simulated 'full stack system' that provides the closest to life experience for training hysteroscopic skills. The simulator is supplied with a standard 4mm hysteroscope, two hysteroscopic instruments, an HD camera system and 15" monitor and monitor stand. The compact table top design makes it the perfect system for use in simulation centres, doctors' offices or as a mobile training solution.

The camera and hysteroscope are powered by the Bozzini light system. The Bozzini light system allows simulated endoscopic procedures to be performed with real endoscopic devices at an affordable price.


Included in the package:

- Hysteroscopy base unit to position uterus
- Clamp for fixing base unit
- 2x non pathological uterus – also used for interchangeable polyps and adhesions
- 1x uterus with submucosal fibroid
- 50x polyps
- 3 rolls of adhesions
- 1x 4mm, 30° or 0° hysteroscope
- 1x hysteroscope introducer
- 1x hysteroscopic grasping forceps
- 1x hysteroscopic scissors
- 1x camera unit
- 1x Bozzini light source
- 1x light lead
- 1x 15" monitor
- 1x monitor stand

Product code: INO-BOZ-HYST

2 Year Extended Warranty


Product code: INO-BOZ-WARR



Key components of the Bozzini Hysteroscopy Simulator

- Uterus base unit and stand

- Bozzini light source




- Hysteroscopic grasping forceps and scissors

- Practise scope handling with real scope and camera



Internal views of pathologies



- 2 x Non pathological uterus and 1 x uterus with submucosal fibroid, polyps and adhesions

Close to life pathologies

Resection of polyps

Resection of adhesions



e: info@inovus.org
t: +44 (0)1744 752 952
w: www.inovus.org



All products contained in this catalogue are Copyright Inovus Ltd 2022
Full price list available on request: info@inovus.org