

# LapAR™ Take Home Simulator

Specification requirements for Laptop or PC to run the LapAR™ take home software.

	<b>Minimum specification</b>	<b>Recommended specification</b>
Processor	Intel I5 CPU	Intel I7 8th Generation CPU
Graphics Card	Nvidia 1050 GPU	Nvidia 1650 GPU
Memory	6GB RAM	16GB RAM
Storage	128GB HDD/SSD	128GB HDD OR SSD

a c c e s s i b l e   a f f o r d a b l e   f u n c t i o n a l

E [info@inovus.org](mailto:info@inovus.org)

T +44 (0)1744 752 952

W [www.inovus.org](http://www.inovus.org)

   @inovusmedical