



Windows 8/10 Camera setup guide:

Step 1: Type 'camera' into the start menu text field and select the camera app.



Step 2: If you can see the inside of the simulator at this stage skip to step 4, if not you will need to switch the computers default camera.



a c c e s s i b l e a f f o r d a b l e f u n c t i o n a l

E info@inovus.org


T +44 (0)1744 752 952

W www.inovus.org


All products are Copyright Inovus Ltd 2020

   @inovusmedical

Step 3: You can do this by clicking the icon circled in red.



Step 4: Operate



a c c e s s i b l e a f f o r d a b l e f u n c t i o n a l

E info@inovus.org

T +44 (0)1744 752 952

W www.inovus.org

All products are Copyright Inovus Ltd 2020

   @inovusmedical